



LOQTEQ® Klavikulaplatten 2.7/3.5 • Clavicle Plates 2.7/3.5

IC 6934-00 Implantateset • Set of Implants



ARTIKEL	ARTICLE	ANZAHL QUANTITY	ART.-NR. ART.-NO.
1 Leersieb Implantate LOQTEQ® Klavikula 2.7/3.5	Tray for implants LOQTEQ® Clavicle 2.7/3.5, empty	1	IC 6934-01
2 Schraubenbank leer, Erweiterungsset LOQTEQ® Klavikula	Screw rack LOQTEQ® Clavicle for extension screw set, empty	1	IC 6934-02
3 Deckel für Leersiebe, groß	Lid for trays, large	1	IC 2008-00
4 Messinstrument für Schrauben 2.5/2.7	Depth gauge for screws 2.5/2.7	1	IS 7903-00
5 Elevatorium schmal, gebogen	Elevatorium small, bended	1	IU 6010-00
6 Bohrer Schnellkupplung Ø2.0, L 110, Wendel 25	Twist drill Ø2.0, L 110, coil 25, quick coupling	2	IU 7420-10
7 Handgriff rund für Schnellkuppl. mit Drehmomentbegr. 1.5 Nm	Handle round with quick coupling, with torque limiter 1.5 Nm	1	IU 7707-00
8 Schraubendrehereinsatz Duo, T8, Schnellkupplung	Screwdriver Duo, T8, quick coupling	1	IU 7815-56
9 Doppelbohrbuchse, Bohrer Ø2.0/2.5	Double drill guide Ø2.0/2.5	1	IU 8125-00
10 Bohrbuchse LOQTEQ® Klavikulaplatten 2.7, hellblau/gold	Drill guide LOQTEQ® Clavicle plates 2.7, light blue/gold	2	IU 8168-20
11 LOQTEQ® AcroPlate® 3.5, 5-Loch, R, Titan	LOQTEQ® AcroPlate 3.5, 5 holes, R, Titanium	1	PK 3511-05-2
12 LOQTEQ® AcroPlate® 3.5, 7-Loch, R, Titan	LOQTEQ® AcroPlate 3.5, 7 holes, R, Titanium	1	PK 3511-07-2
13 LOQTEQ® AcroPlate® 3.5, 5-Loch, L, Titan	LOQTEQ® AcroPlate 3.5, 5 holes, L, Titanium	1	PK 3512-05-2
14 LOQTEQ® AcroPlate® 3.5, 7-Loch, L, Titan	LOQTEQ® AcroPlate 3.5, 7 holes, L, Titanium	1	PK 3512-07-2
15 LOQTEQ® Klavikulaschaftplatte 3.5, 6-Loch, L 76, R, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 6 holes, L 76, R, Titanium	1	PK 3521-06-2
16 LOQTEQ® Klavikulaschaftplatte 3.5, 7-Loch, L 88, R, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 7 holes, L 88, R, Titanium	1	PK 3521-07-2
17 LOQTEQ® Klavikulaschaftplatte 3.5, 8-Loch, L 95, R, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 8 holes, L 95, R, Titanium	1	PK 3521-08-2
18 LOQTEQ® Klavikulaschaftplatte 3.5, 10-Loch, L 115, R, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 10 holes, L 115, R, Titanium	1	PK 3521-10-2
19 LOQTEQ® Klavikulaschaftplatte 3.5, 6-Loch, L 76, L, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 6 holes, L 76, L, Titanium	1	PK 3522-06-2
20 LOQTEQ® Klavikulaschaftplatte 3.5, 7-Loch, L 88, L, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 7 holes, L 88, L, Titanium	1	PK 3522-07-2
21 LOQTEQ® Klavikulaschaftplatte 3.5, 8-Loch, L 95, L, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 8 holes, L 95, L, Titanium	1	PK 3522-08-2
22 LOQTEQ® Klavikulaschaftplatte 3.5, 10-Loch, L 115, L, Titan	LOQTEQ® Clavicle Shaft Plate 3.5, 10 holes, L 115, L, Titanium	1	PK 3522-10-2
23 LOQTEQ® Sup. Lat. Klavikulaplatte 2.7/3.5, 6/4-Loch, L 63, R, Titan	LOQTEQ® Sup. Lat. Clavicle Plate 2.7/3.5, 6/4 holes, L 63, R, Titanium	1	PK 3531-04-2
24 LOQTEQ® Sup. Lat. Klavikulaplatte 2.7/3.5, 6/7-Loch, L 93, R, Titan	LOQTEQ® Sup. Lat. Clavicle Plate 2.7/3.5, 6/7 holes, L 93, R, Titanium	1	PK 3531-07-2
25 LOQTEQ® Sup. Lat. Klavikulaplatte 2.7/3.5, 6/4-Loch, L 63, L, Titan	LOQTEQ® Sup. Lat. Clavicle Plate 2.7/3.5, 6/4 holes, L 63, L, Titanium	1	PK 3532-04-2
26 LOQTEQ® Sup. Lat. Klavikulaplatte 2.7/3.5, 6/7-Loch, L 93, L, Titan	LOQTEQ® Sup. Lat. Clavicle Plate 2.7/3.5, 6/7 holes, L 93, R, Titanium	1	PK 3532-07-2
27 Kortikalisschraube 2.5, kl. Kopf T8, selbstschn. L 10-28, Titan	Cortical Screw 2.5, small head T8, self-tapp. L 10-28, Titanium	0	SK 2512-xx-2
28 LOQTEQ® Kortikalisschr. 2.7, kl. Kopf T8, selbstschn. L 10-26, Titan	LOQTEQ® Cortical Screw 2.7, small head T8, self-tapp. L 10-26, Titanium	5	SK 2726-xx-2